# Pelvic Fixation of Growing Rods



Paul Sponseller MD, Behrooz Akbarnia MD, George Thompson MD, Richard McCarthy MD, John Emans MD, Marc Asher MD, Muharrem Yazici MD, David Skaggs MD., Connie Poe-Kochert CRNP, Pat Kostial RN, BSN, Tina, Chen, BS

- for the GSSG

#### **Financial Disclosures**

#### Presenter: Paul Sponseller - (a,b) DePuy Spine

- Co-Authors:
- George Thompson
- Behrooz Akbarnia (a) DePuy Spine
- John Emans
- Marc Asher

- a. Grants/Research Support
- b. Consultant
- c. Stock/Shareholder
- d. Speakers' Bureau
- e. Other Financial Support

#### Introduction

- Growing rods = "internal brace" for young patients with severe scoliosis extending to the pelvis
- No prior studies of growing rods to pelvis
  - How do foundations behave over time?
- This project analyzed the outcomes and complications unique to this construct

#### **Patients and Methods**

- 22 patients from 8 centers
- Indications/Inclusion criteria
  - Severe pelvic obliquity
  - Distal deformity
  - Lack of satisfactory alternative anchor sites
  - ≥ 2 years treatment with growing rods fixed to the pelvis

## Diagnoses

- Myelomeningocele 4
- Congenital 4
- Cerebral palsy 3
- Arthrogryposis 1
- SMA 1
- Miscellaneous/syndromic 9

#### **Patient Characteristics**

- Age at surgery 6.1 + 3.1 years.
- Preop curve 86 ± 22°
- Coronal imbalance 9.7 ± 8.2 cm
- Follow up 50 months.
- Mean of 2.9 + 1.8 lengthenings (1-9)

### Anchor types

- Pelvic Fixation:
  - Iliac screws or rods -17
    - 3 s-rods
    - 5 iliac rods
    - 11 screws
  - Sacral hooks -5
- Proximal fixation
  - Hooks-12
  - Screws-10
- Dual rods used in 18 patients; single in 4
- Distal crosslink used to improve stability

#### Results

 Curve improved from 86° to 47° ±19° at final follow up.

- Correction same w. Sacral vs iliac anchors
  - -45% vs 46%; ns
- Coronal imbalance improved from 9.7 to 4.4 cm

#### Results

- Mean increase in T1-S1 length 7.3 ± 2.9
  cm during distraction
- Seven patients have undergone final fusion at a mean of 10.8 + 1.4 years

# CP 4 yrs; final fusion







Arthrogryposis




Pre-op →





1st Distraction → Last Follow-up



Syndromic Curve








# SMA II









### Complications

- 6 deep wound infections
- 6 distal fixation complications; all salvaged
- 3 rod breakages
  - this rate did not differ statistically from the rate for dual growing rods as a whole
  - (3/22 vs/ 25/144; ns)

# Complications by distal anchor type

- Hooks:
  - 2/3 distal erosions
- S-rods
  - 2/3 migrations
- Screws:
  - 2/11 edge prominences
- Iliac rods
  - no complications reported

## Suggestions

- Bilateral pelvic fixation
- Maximize intra-pelvic length
- Rigid distal cross link
- Avoid prominent implant edges

#### Conclusions

- Pelvic fixation an effective caudal foundation for growing rods
- Both screws and hooks satisfactory;
  - Iliac rods/screws had lower complication rate
- Age & Indications still being defined
  - For start and stop
- 10+ cm realistic
- Indicated when definitive fusion will involve pelvis

## Thank You

