Some Connectors in Distraction Based Growing Rods Fail More Than Others

Christopher Lee, BS; Karen S. Myung, MD, PhD; David L. Skaggs, MD

ICEOS 2011

Correspondence: kmyung@chla.usc.edu



Dislosures

- Christopher Lee None
- Karen Myung None
- David Skaggs Medtronic (b,d); Stryker (b,d); AO
 North America (d); Biomet (e)



Hypothesis

• Some connectors fail more than others.



better than







Background

- Distraction-based growing rods have ~2 complications per patient¹
- No study has examined connector failures in growing rod constructs.



Methods

- 30 consecutive patients identified retrospectively at a single institution.
- All patients had distraction-based growing rods.
- Only connectors subject to distractive forces included.
- All patients had minimum 2 year followup (range, 24-83 months).



Connector Failures

- Total of 12% failure =12/98 connectors:
 - All connectors subject to distractive forces.

 Mean time to failure = after 3rd lengthening





Connector Failures

• Dual-rod connectors vs. single-rod connectors:

- Dual-rod connectors: 32% failure (6/19 constructs)
- Single-rod connectors: 43% failure (6/14 constructs)
- No difference in failure rates (p=0.38).



2 Types of Interface with Rods

<u>Circular Lock</u> Older design 41% failure rate V-shape Lock Previously patent protected 4% failure rate





NIVERSIT

Type of Connector	Failure Rate (n, %)	Connectors
<u>Circular Lock</u> - Medtronic® Wedding Band	7/17, 41%	
<u>V shaped Lock -</u> Medtronic [®] and Isola [®]	1/27, 4%	
Stryker [®] , Open	3/8, 38%	
Longitudinal	1/34, 3%	

28% less T1-S1 Growth with connector failure

- Connector Failure = 6.5 mm/yr
- No failure = 9.0 mm/yr
- (p>0.05) but underpowered
- Seems logical that if distraction is not present, spine growth not driven



Discussion

- Some connectors fail more than others:
 - side-to-side connectors with tapered Vshaped locks failed 10 times less often than those with circular locks.
 - "Wedding Bands" with circular holes fail almost as much as first marriages (41% vs. >50%)
- Age of initial instrumentation may matter:
 - >5 years old: 0 connector failures
 - Greater complications reported in patients treated at younger age^{2,3}



Take Home Lesson: Uselongitudinal connectors or sideto-side connectors with V- locks









CHILDREN'S ORTHOPAEDIC CENTER