

SHILLA: Techniques, Tips and Tricks

Scott J. Luhmann, M.D.

Associate Professor

Department of Orthopaedic Surgery

Washington University School of Medicine

Chief of Spine Surgery, Shriner's Hospital

St. Louis, Missouri

7th ICEOS, November 22, 2013



Washington University in St. Louis

SCHOOL OF MEDICINE




Disclosures


- Speaker's bureau for:
 - Medtronic Sofamor Danek
 - Stryker Spine
- Consultant for:
 - Medtronic Sofamor Danek
 - Stryker Spine
 - Orthofix
 - Depuy Synthes
- Royalties:
 - Globus Medical

SHILLA

- 3 topics
 - Metallosis
 - Rod breakage
 - Implant prominence



Pictures: Rick McCarthy




Metallolysis:

Shilla growing rods in a caprine animal model: a pilot study

McCarthy RE, Sucato D, Turner JL , et al
CORR 2010

5.5 SS




- 11 2 m/o immature goats
- Spines harvested at 6 m p-op
- Mean axial growth over construct: 48mm
- No apical spinal stenosis

Shilla growing rods in a caprine animal model: a pilot study

McCarthy RE, Sucato D, Turner JL , et al
CORR 2010



- No implant failure
- Minor wear at rod/screw



5.5 SS



2+9 years post-implantation

Wear:
Dorsal > Ventral




Wear:
Caudal > Cephalad

Wear Patterns:
Rods matched set screws




2+9 years post-implantation




Shilla growing rods in a caprine animal model: a pilot study

McCarthy RE, Sucato D, Turner JL , et al
CORR 2010

- Posterior
 - SHILLA screws (cephalad and caudal)
 - Metallic debris in soft tissues
 - Moderate to extensive inflammatory reaction
 - Apical (fused) segment: no debris



Intraoperative




Pictures: Rick McCarthy

Shilla growing rods in a caprine animal model: a pilot study

McCarthy RE, Sucato D, Turner JL , et al
CORR 2010

- Anterior lymph nodes
 - SHILLA screws (caudal): Metallic debris
 - Apical (fused) segment: no debris
- Systemic dissemination: unknown
- Not a synovial joint
- Design improvements



SHILLA: Rod Breakage


- Do Shilla Rods Have Acceptable Complications and Fewer Surgeries?

McCarthy RE, Luhmann SJ, Lenke LG, McCullough FL;
SRS 2009


- **Rod breakage more common in:**

- Active, ambulatory children
- 3.5 mm SS rod

7 y/o
3.5 mm




12 y/o
4.5 mm




- 4.5 better than 3.5
- Is 5.5 better than 4.5?
- Would rather replace broken rods than failure of bone-screw interface.

SHILLA: Rod Breakage


14.5 yrs
4.5 yrs post-SHILLA



7 y/o




12 y/o




Majority of rods break caudal to apical fusion

Apical Fusion


2 level





3 levels



4 levels



Intercalary SHILLA Screws




SHILLA: Rod Breakage

7 y/o male
NF

3-level apical
fusion

6 months postop




SHILLA: Implant Prominence

- Fully seat all pedicle screws
- Larger rod size means more prominence
- Rod contouring:
 - Cephalad: bend into more kyphosis
 - Caudal: bend into lordosis
- Leave 1 inch of rod at top and bottom of construct.




5 ½ years s/p T3-L3 Shilla: 1 reoperation


6 y/o male
6 wks
postop




3 y postop



4 ½ y postop



s/p revision
11+5 yrs



Implant Prominence

5.5 SHILLA

4.5 SHILLA

5.5 CMAS



17mm

12mm



16mm

12mm



15mm

10mm

Thank you

