#### **ICEOS 2016**

#### Five or more Proximal Anchors and including the Upper End Vertebrae (UEV) Protects Against Reoperation in Growth Friendly Constructs

Liam Harris BS<sup>1</sup>, Lindsay M Andras MD<sup>1</sup>, Gregory M Mundis MD<sup>2</sup>, Paul D Sponseller MD MBA<sup>3</sup>, John B Emans MD<sup>4</sup>, David L Skaggs MD MMM<sup>1</sup>

1.Children's Orthopaedic Center, Children's Hospital Los Angeles; 2. San Diego Center for Spinal Disorders 3. Department of Orthopaedic Surgery, Johns Hopkins Childrens Hospital, Johns Hopkins University; 4.Department of Orthopaedic Surgery, Boston Children's Hospital, Harvard Medical School;

USC University of Southern California



#### Disclosures

- a. Grants/Research Support
- b. Consultant
- c. Stock/Shareholder
- d. Speakers' Bureau
- e. Editorial/Governing Board
- f. Other Financial Support

- Liam Harris BS None
- Lindsay M. Andras MD Eli Lilly (c); Orthobullets (f)
- Paul D. Sponseller MD, MBA -DePuy (a, b, f); Globus Medical (f); Journal of Bone and Joint Surgery (e); Journal of Bone and Joint Surgery Oakstone Medical (f); Scoliosis Research Society (e)
- John Emans MD Journal of Children's Orthopedics (e); Medtronic (b); Synthes (b, f)
- David L. Skaggs MD MMM SRS & POSNA (a); Biomet; Medtronic; Orthobullets; Grand Rounds (b); Zipline Medical, Inc. (b & c); Biomet; Medtronic; Johnson & Johnson (d); Wolters Kluwer Health - Lippincott Williams & Wilkins; Biomet Spine (f)





# Background

- Implant failure in Growing Spine constructs can approach 80%<sup>3</sup>
- Common reasons for proximal revision:
  - Proximal Junctional Kyphosis<sup>4</sup>
  - Proximal Anchor Pullout<sup>5</sup>
- Biomechanical studies suggest proximal anchor type, number and location may influence implant failure rates <sup>6</sup>







# Purpose

To assess the impact of number, type and placement of proximal anchors on complication and revision rates following surgery with growth friendly constructs



USC University of Southern California



CHILDREN'S DRTHOPAEDIC CENTER

## Methods

PA

- Retrospective review, multicenter database
- Index instrumentation at < 10 years of age with <u>></u> 2 year follow up
- Upper end vertebrae of kyphosis (UEV) recorded on lateral radiographs
- Statistical significance analyzed via t-test and multilinear regression analysis







#### Results

# 357 patients mean f/u of 6 years: -306 growing rods with spine anchors -16 growing rods with rib anchors

- -35 VEPTRs
- Implant type (VEPTR vs. Growing Rod) was not associated with Anchor Pullout (p=0.150)







## **Results- Anchor Pullout**

- 22% (77/357) patients anchor pullout
- Type of Anchor (screw vs spine hook vs rib hook) did not affect pullout rate (p=0.850)

Anchors	Pullout
<u>&gt;</u> 2	22%
<u>&gt;</u> 3	20%
<u>&gt;</u> 4	20%
<u>&gt;</u> 5	12%
<u>&gt;</u> 6	10%



# 2,3,or 4 Vs. 5 or 6 anchors p=0.01

#### **Results - Proximal Revision**

 Mean Index Proximal Instrumentation : T2 (range C7- T9)

 Instrumentation below UEV was associated with increased proximal revision (p=0.026)



**UEV** of

Sagittal

#### **Results - Proximal Extension**

- 24% (47/200) Proximal Extension of Implants (final fusion or revision growth friendly implants)
  - Index Instrumented at or above UEV: **20%** (26/130)
  - Index Instrumented below UEV: 30%(21/70)
  - Index Instrumentation below UEV associated with higher rate of proximal revision (p=0.027)





#### **Results - Proximal Revision**

#### • Not associated with Proximal Revision:

- Cobb Angle
- Total Anchors
- Type of Anchor
- Type of Implant
- Number of vertebrae instrumented
- Level of UEV of Kyphosis





CENTER

#### Conclusions

- 1.  $\geq$ 5 anchors less pullout (p=0.010)
- 2. Anchors at or above UEV (Sagittal) less proximal extension (p=0.026)

- 3. Not Significant:
  - Anchor Type (spine hook, rib hook, screw)
     Implant Type (MEDTD or CD)
  - 2. Implant Type (VEPTR or GR)



Southern California



#### References

- 1. McCarthy, R.E., et al., *The Shilla growth guidance technique for early-onset spinal deformities at 2-year follow-up: a preliminary report*. J Pediatr Orthop, 2014. **34**(1): p. 1-7.
- 2. La Rosa, G., L. Oggiano, and L. Ruzzini, *Magnetically Controlled* Growing Rods for the Management of Early-onset Scoliosis: A Preliminary Report. J Pediatr Orthop, 2015
- 3. Sankar, W.N., D.C. Acevedo, and D.L. Skaggs, Comparison of complications among growing spinal implants. Spine (Phila Pa 1976), 2010. 35(23): p. 2091-6.
- 4. El-Hawary, R., et al., What is the Risk of Developing Proximal Junctional Kyphosis During Growth Friendly Treatments for Early-onset Scoliosis? J Pediatr Orthop, 2015.
- 5. Watanabe, K., et al., Risk factors for complications associated with growing-rod surgery for early-onset scoliosis. Spine (Phila Pa 1976), 2013. 38(8): p. E464-8.
- 6. Mahar, A.T., et al., Biomechanical comparison of different anchors (foundations) for the pediatric dual growing rod technique. Spine J, 2008. 8(6): p. 933-9.





CHILDREN'S

•11